

Introduction

SBB employs an innovative sequencing chemistry that delivers extraordinary accuracy due to ultra-low levels of sequencing noise. By incorporating native nucleotides, optimizing chemistry at each step, and minimizing molecular scarring, SBB provides significantly higher quality data than existing technologies.

Sequencing by binding (SBB)



- 2 Interrogate TACGAGT
 ADGODOAGDE
- 3 Activate TACGAG \rightarrow ADGGTGAGDT
- 4 Incorporate TAGGAGT

Figure 1. Conceptual diagram of SBB sequencing

Key features of SBB

- Extraordinary accuracy 90% of bases ≥Q40
- · Ultra-low sequencing noise
- High-confidence sequencing through repetitive regions
- · Negligible duplication rate
- Minimal index-hopping

(1) Each cycle initiates with a 3' reversible blocked nucleotide. (2) Fluorescently labelled nucleotides are then flowed over the flow cell allowing for the appropriate base to bind. Unbound nucleotides are washed away so that the base can be interrogated with reduced background signal. (3) The 3' end of the nucleotide is activated via removal of the reversible terminator. (4) Native, unlabeled, reversible blocked nucleotides are flowed over the flow cell and the cognate base can incorporate into the growing strand. The process is then repeated for each new sequenced base.

By separating the interrogation and incorporation steps, and optimizing chemistries for each step, SBB chemistry enables high signal-to-noise levels to be achieved. Moreover, the lack of molecular scarring from the use of native nucleotides for incorporation drives industry-leading accuracy.

	Sample applications		
	WGS ¹ (2x150bp)	scRNA-Seq ² (28×90bp)	Exome ² (2x100bp)
Number of paired-end reads	1,033 M	1,000 M	1,062 M
Total yield	151 Gb	58 Gb	116 Gb
Q score (≥90% of bases)	50	46	49

^{1.} Sequenced with Onso native library prep kit

Table 1. Sample output metrics from the Onso[™] sequencing system.




Figure 2. SBB sequencing delivers near perfect reads with >90% of bases above Q40, low duplication rates, minimal index-hopping, and sequencing through difficult and repetitive regions.




Figure 3. The Onso sequencing system produces near perfect reads through the variable poly-T region within the *TOMM40* gene. (a) SBS sequencing was unable to fully resolve the low-complexity region, with many false positive calls. (b) SBB sequencing was able to fully sequence through and confidently resolve this region, resulting in a call of T16 vs the original estimate of T35.

^{2.} Sequenced using converted libraries

SBB complete workflow



Onso library prep

- Employ standardized NGS library prep optimized for SBB sequencing
- Create Onso libraries directly from samples
- Adapt existing P5/P7 libraries to an Onso-compatible format



Onso sequencing

- Easily implement workflow into new and existing NGS labs
- Produce highly accurate reads with SBB chemistry
- Individually addressable flow cell lanes



Data analysis

- Generate industry-standard FASTQ output
- Easily integrate with currently available secondary analysis tools



Explore the applications: pacb.com/applications



Connect with PacBio for more info: North America: nasales@pacb.com South America: sasales@pacb.com

EMEA: emea@pacb.com

Asia Pacific: apsales@pacb.com



Contact a certified service provider pacb.com/CSP

Research use only. Not for use in diagnostic procedures. © 2024 Pacific Biosciences of California, Inc. ("PacBio"). All rights reserved. Information in this document is subject to change without notice. PacBio assumes no responsibility for any errors or omissions in this document. Certain notices, terms, conditions and/or use restrictions may pertain to your use of PacBio products and/or third-party products. Refer to the applicable PacBio terms and conditions of sale and to the applicable license terms at pacb.com/license. Pacific Biosciences, the PacBio logo, PacBio, Circulomics, Omniome, SMRT, SMRTbell, Iso-Seq, Sequel, Nanobind, SBB, Revio, Onso, Apton, Kinnex, and PureTarget are trademarks of PacBio.

