


Iso-Seq on the Sequel II System

- Generate full-length transcripts of 10 kb or longer
- High accuracy (>99%) for ORF prediction
- No reference genome required
- Bioinformatics tools from raw data to functional annotation
- Many applications, including: [1]
 - Genome annotation
 - Novel gene and isoform discovery
 - Fusion gene detection
 - Allele-specific isoform expression analysis
 - Better reference for RNA-seq quantification
 - Assess genome assembly quality
 - Single cell analysis

SQANTI2 for Complex Transcript Classification




SQANTI2 classifies isoforms from an Alzheimer brain Iso-Seq dataset at a complex locus.

Novel isoforms (NIC) and readthrough transcripts of multiple genes identified.

Iso-Seq Express Kit [2]

- Input 60-300 ng total RNA
- Full-length cDNA
- Multiplexing support




Sequel II System

- 1 SMRT Cell 8M for whole transcriptome
- Up to 4 million full-length reads

Main Bioinformatics Tools

	Input	Output
Sequencing		subreads .bam
Iso-Seq Analysis	subreads .bam or ccs .bam	Collapsed unique transcripts (GFF, FASTA)
Transcript Classification	Unique transcripts Reference genome Annotation (GTF) CAGE Peak Junction data...	Transcript classification Junction classification Figures
Functional Annotation	SQANTI output	Annotated GTF
Differential Analysis	Experimental design Annotated GTFs	

SQANTI2 Transcript Classification



Supporting Bioinformatics Tools

	<ul style="list-style-type: none"> - collapse redundant transcripts - merge multi-sample output - saturation curve - file format conversion - single cell analysis 		<ul style="list-style-type: none"> - gene family finding - genome reconstruction - evaluate assembly
	<ul style="list-style-type: none"> - collapse redundant transcripts - merge multi-sample output - NMD/ORF prediction - transcript filtering 		<ul style="list-style-type: none"> - long read processing & annotation pipeline developed independently by ENCODE4

	FL Reads	Unique Genes	Unique Transcripts
UHRR [3]	4,734,362	16,328	183,689
Alzheimer Brain [4]	4,277,293	17,670	162,290

[1] Low et al., "Haplotype-Resolved Cattle Genomes Provide Insights Into Structural Variation and Adaptation", *bioRxiv* (2019)
 Beiki et al., "Improved annotation of the domestic pig genome through integration of Iso-Seq and RNA-seq data", *BMC Genomics* (2019)
 Wang et al., "Reviving the Transcriptome Studies: An Insight Into the Emergence of Single-Molecule Transcriptome Sequencing", *Front Genet* (2019)
 [2] Iso-Seq Express: <https://www.pacb.com/wp-content/uploads/Procedure-Checklist-Iso-Seq-Express-Template-Preparation-for-Sequel-II-Systems.pdf>
 [3] UHRR: [https://github.com/PacificBiosciences/DevNet/wiki/Sequel-II-System-Data-Release:-Universal-Human-Reference-\(UHR\)-Iso-Seq](https://github.com/PacificBiosciences/DevNet/wiki/Sequel-II-System-Data-Release:-Universal-Human-Reference-(UHR)-Iso-Seq)
 [4] Alzheimer brain: https://downloads.pacbloud.com/public/dataset/Alzheimer2019_IsoSeq/

